

GS



Growth and versatility that keeps up with busy kids. Named for its built-in growth system, the Zippie® GS™ offers frame and component growth, standard and reverse wheel configurations, and folding and rigid frame options to keep up with children's ever-changing lives!

Growth capacity with 3" depth built into the frame. **Image 1.**

Ultra Light Aerospace 7000 Aluminum Frame. **Image 2.**

Frame, fork and deck, axle, wheel and guard color options. **Image 3.**

A custom Ride Designs® seating system. **Image 4.**



Image 1




Image 2



Image 3




Image 4

Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.


SPECIFICATIONS

ITEM	EIZGS
Frame	Aluminum
Frame type	Folding / Rigid
Seat Width	10" a 18" / 25.4 cm a 45.7 cm
Seat Depth	10" a 20" / 25.4 cm a 50.8 cm
User Weight Capacity	75 kg / 165 lbs
Product Weight	13.2 kg / 29 lbs
Accessories	Optional
Settings	In measurements

Lightweight and foldable frame easy to transport and store



Available colours:



BENEFITS

Technology of a rigid frame and folding frame for active and growing children.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com



ver_2025.11.11E

