

# GottaGo Portable Toilet Seat



Image 1



Image 2



Image 3



Image 4

Innovative portable toilet that facilitates toilet training for families of children with special needs.

Portable and height adjustable, making it suitable for any environment at home or away from home. **Image 1.**

Fastening straps for different parts or user needs, practical and safe. **Image 2.**

Backpack capable of holding your GottaGo and everything that comes with it, it also has space for wipes, diapers, spare clothes. **Image 3.**

Potty insert for use when toilets are not available, easily attaches to the GottaGo seat and comes with its own washable carrying bag. **Image 4.**



# Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

## SPECIFICATIONS

Sizes	1 y 2
Age approx.	(1) 2 to 5 years   (2) 4 to 9 years
Hip width	(1) 22.5 cm   (2) 27 cm
Deep seat	(1) 22 cm   (2) 29 cm
Backrest inclination	(1) 8°   (2) 25°
Thigh angle	(1) 40°   (2) 40°
Backrest height	(1) 24 cm   (2) 36 cm
Max user weight	(1) 30 kg   (2) 30 kg
Product weight	(1) 5.8 kg   (2) 6.1 kg
Seat height	(1) 30 cm   (2) 46 cm

2 sizes available



## BENEFITS

Suitable for patients with neuronal and psychomotor damage or cognitive psychic disability.

## Discover our products



### DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

