

Hydraulic Patient Lift



The Hydraulic Patient Lift from ProBasics is economically priced and has a 450 lb. weight capacity. The dual purpose padded swivel bar can use 4 or 6 point slings. The rotating handle has a manual hydraulic pump and the frame is a chip resistant brown vein finish.

Heavy-duty round tube steel construction with a scratch-resistant brown finish. **Image 1.**

Made with water-resistant materials and easy to clean. **Image 2.**

Heavy duty 5" solid wheels (2) with reinforced locking. **Image 3.**

Various sizes of nylon four-point strap washable and can be used in wet or dry environments. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS	
ITEM	PBPLSM001
Closed base	56 cm / 22"
Open base	108 cm / 42.25"
Maximum height	162 cm / 64"
Minimum height	92 cm / 22"
Crane height	198 cm / 78"
Wheel size	5"
Weight capacity	204 kg / 450 lbs
Product weight	31 kg / 68 lbs
System	Hydraulic manual
Accessories	Transfer meshes

BENEFITS

Essential in transfers for patients and caregivers in a comfortable and safe way.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

