

Kid Kart Xpress



zipie®

The KID-KART® Xpress early intervention mobility kids wheelchair features superior positioning, contemporary style, and user-friendly functionality for children 0-3 years old weighing up to 55 lbs.

Removable Seat and Folding Base.
Image 1.

With 20° of forward tilt and 45° of rearward tilt a 65° adjustable tilt range.
Image 2.

Adjustable Positioning Components.
Image 3.

The Fixed Height Feeder Base and Height Adjustable Booster Base keep your child properly positioned during feeding, activities, or therapy. **Image 4.**



Image 1



Image 2



Image 3



Image 4

Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS

ITEM	EIXP
Seat width	6" to 12.5"
Seat depth	5.5" to 11"
Product weight	15.4 kg / 34 lbs
Backrest height	13" to 22"
Weight capacity	25 kg / 55 lbs
Incline range	65° (20° anterior, 45° posterior)
front wheels	7" and 8"
Rear wheels	10" - 14"
Accessories	Optional
Colors	On base and hoods

Required user
measurement
(Total height)



BENEFITS

A special needs stroller that provides mild to moderate positioning support

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com



ver_2025.11.11E

