

Gait Trainer Klip



It is a lightweight rear gait trainer for children, teenagers and young adults with special needs who require support, assistance in their daily walking. Promotes walking in an upright position.

Folding frame easy to carry and store. **Image 1.**

Adjustable height and folding seat option, to adapt to the growth of the user. **Image 2.**

Wheel lock with an audible "click" for easier control and Foot-activated swivel. **Image 3.**

Adjust the resistance on the rear wheels to control speed. **Image 4.**



Image 1



Image 2



Image 3



Image 4




Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Frame	Blue and red aluminum
Handle height	40.6 - 94cm / 16" - 37"
Width of handles	34.3 - 43.2cm / 13.5" - 17"
Base Width	58.4 - 94cm / 23" - 37"
Weight capacity	86 kg / 190 lbs
Product weight	From 4 kg / From 9 lbs

Required user
measurement
(Elbow to the floor)



BENEFITS

For users who need to improve their ambulatory skills and develop weight support

Discover our products



DESIGNED FOR OUR TERRITORIES



Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

