

MOSAIC Cushion

ROHO 



This innovative design conforms to your body's contours, offering a comfortable sinking sensation. Available in a wide range of sizes to ensure the best fit, it offers basic protection for skin and soft tissues.

The stretch and antibacterial cover adapts and maintains to any surface or wheelchair. **Image 1.**

The unique air cell design has the ability to adjust and adapt to the shape of the user's body. **Image 2.**

Made of high quality and ultralight polyvinyl in various sizes. **Image 3.**

Water resistant and easy to clean materials. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

Wheelchair cushions designed to provide the user with support, comfort and protection.

SPECIFICATIONS	
Width range	16.25" to 20.25" 41.5 cm to 51.5 cm
Length	16.25" to 18.25" 41.5 cm to 46.5 cm
Product weight	0.4 kg / 0.88 lbs
Height	7.6 cm / 3"
Weight capacity	143 kg / 315 lbs
Outer cover	Elastic cover
Air pump	Included (manual)



BENEFITS

Designed for users at risk of skin degradation requiring mild stability

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

