

Mobile Stander

 rifton.



Eye to eye, smile to smile – engaging with peers is so much more rewarding for individuals with special needs when they can interact in an upright and mobile position. The Rifton Mobile Stander, a mobile prone stander, makes it easy.

Versatile arm prompts include height, rotation, angle, lateral proximity and forward/backward adjustments.

Image 1.

The large wheels with brakes are placed with just a button. Spoke guards protect fingers and rims keep hands clean.

The conversion kit converts the smallest to the next size up.

The communication tray is fully adjustable and holds small objects or activities.



Image 1



Image 2



Image 3



Image 4




Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Frame measurements	64 to 165 cm / 25" to 65"
Max weight of user	23 to 79 kg / 50 up to 175 lbs
Product weight	12.7 - 143 kg / 28 - 65 lbs
Accessories	Optional
Wheel diameter	20", 27" and 35"
Settings	No tool required
Frame color	Blue, Lime and Pink

Required user
measurement
(Total height)



BENEFITS

Allows upright position for adults and children with physical disabilities.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

