

# Q200 R



Ultra-compact RWD power wheelchair with super-tight indoor maneuverability and ZERO compromise on outdoor stability and performance.

The Captain Seat provides a relaxed sitting position and features a semi-recline fold down back. **Image 1.**

Available with a 22° tilt and power centermount leg support. **Image 2.**

SMART base technology. **Image 3.**

Get expressive with color! **Image 4.**



Image 1




Image 2



Image 3




Image 4



# Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

Seat Width 15", 16", 18", 20" & 21"

Seat Depth 16" to 22"

Backrest Height 18" to 30"

User Weight capacity 136 kg / 300 lbs

Product weight From 126 kg / 278 lbs

Turning radius 18.3" without footrest

Front wheels 9"

Rear wheels 12" or 14"

Speed 9.6 kmh / 6 mph

Available colours



Black



Blue



Green



Orange



Red



White



Silver

2 sealed batteries



## BENEFITS

For active people mainly spinal and motor injuries.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

