

Q50 R



It is your perfect travel companion and for daily activities. Equipped with a lithium battery approved for air travel, it offers a range of up to 7.5 miles and allows its range to be doubled up to 15 miles.

Weighing just 32 lbs. (without batteries) and made with top quality carbon fiber. **Image 1.**

The ultralight backrest features breathable upholstery, allowing for optimal airflow and temperature regulation. **Image 2.**

Foldable in seconds, it also stores comfortably in the tightest places without removing the battery. **Image 3.**

Reliable and established nVR2 brushless motor controller responds precisely to your commands and battery charging. **Image 4.**




Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	17"
Seat depth	16.5"
Backrest height	18"
Product weight	(without battery) 14.5 kg 32 lbs
Max user weight	135 kg 300 lbs
Rear wheels	8"
Caster	6.5"
Maximum speed	6 kmh 3.7 mph
Battery	24V, 10Ah Lithium
Motor	High power 250 W, brushless
Overall dimensions	(L/W/H) 36.5"/22.4"/36.3"
Folded dimensions	(L/W/H) 30.2"/22.4"/14.6"

Charge
Even When Folded



Superbright
running light

5 eye popping colors



BENEFITS

For users with partial or total reduced mobility of lower limbs.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES




Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

