

QS5 X



The QUICKIE QS5 X, redefines easy folding thanks to FreeFold - an innovative cross-brace design that reduces folding forces and saves you over 40% energy with every fold! The QS5 X makes light work of lifting, with an (optional) folding fastener that keeps the wheelchair tightly locked to prevent any unintentional opening during transport.

The QS5 X is easy to fold, and it stays folded - thanks to the innovative FreeFold feature. **Image 1.**

Swing-Away Folding chair frame front option. **Image 2.**

Fixed chair frame front option. **Image 3.**

The QS5 X uses a press fit to connect the frame to the oval tube reinforced cross member and cross brace. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Frame	7000 Aluminum Folding
Seat Width	13" to 20"
Seat Depth	13" to 20"
Product Weight	From 7.4 kg / 16.2 lbs
Weight Capacity	135 kg / 300 lbs
Frame Angle	70°, 75°, 80°, 85° and 90°
Center of Gravity	Adjustable 1" - 3"
Accessories	Optional
Adjustments	In sizes
Available Colors	32 available



BENEFITS

For active people mainly spinal and motor injuries.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

