

# Quickie IRIS



The best space tilt wheelchair, the wheelchair with manual tilt up to 55°, the widest range of positioning options, numerous adjustments, and is one of the lightest space wheelchairs available today.

It offers the recommended option for JAY cushions and backrests tailored to the user's specific needs. **Image 1.**

The power tilt option gives the user more control over their body position in space. **Image 2.**

Adjustable and removable frame and accessories for easy transport. **Image 3.**

Features 40° or 55° rotation ranges that can be configured to achieve up to 45° or 60°. **Image 4.**



Image 1



Image 2



Image 3




Image 4



# Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

## SPECIFICATIONS

ITEM	EIZ4-2
Frame	Aluminum
Seat width	14" to 22"
Seat depth	14" to 22"
Weight capacity	135 kg / 300 lbs
Tilting degrees	-5° up to 55°
Tilting system	Tilt in space
Accessories	Optional
Settings	In measurements
Available colours:	

Pedal system to activate tilt and brake	
JAY Lower Extremity Support	
Intuitive foot support designs	
The MONO Single Post Backup System	

## BENEFITS

A special needs chair that provides tilting and positioning support.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

