

Scoot



It is designed to help children with special needs play, explore and participate in a way that best suits their physical and cognitive abilities.

Scoot Backrest provides more support and reduces the amount of trunk control needed. **Image 1.**

Push Pole. **Image 2.**

The footplate can be easily adjusted to your child's exact leg length. **Image 3.**

Advanced Backrest Cushion. This optional extra that gives a supportive 'hug' to your child's trunk. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

Equipment manufactured and adapted for children and adults with disabilities.

SPECIFICATIONS

Backrest height	26 cm / 10.5"
Seat width	23 cm / 9"
Deep seat	40 cm / 15.5"
Backrest width	26 cm / 10.5"
Max user weight	22 kg / 48.5 lbs
Age range	2 to 6 years
Wheel size	9"
Pole height	80 - 89 cm / 31.5" - 35"



BENEFITS

Suitable for children with neuronal and psychomotor damage or cognitive psychic disability.

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

