

# Soft-Touch Sitter



The Sitter is an adaptive seating system that provides support in a variety of sitting positions. Comes in three different variations: Sitter only, fixed wedge base or base with wheels.

The Sitter has anti-microbial properties, is impermeable to fluids, as well as latex-free and peel- and tear-resistant.  
**Image 1.**

Floor Wedge. **Image 2.**

Stationary Base. **Image 3.**

Mobile Base. **Image 4.**



Image 1



Image 2



Image 3








Image 4



# Keeping you in mind. More options.

Clinically sound, innovative design and exceptional value.

SPECIFICATIONS	
User height	Size 1 to 5 75 - 185cm / 30 - 74"
User weight	Size 1 to 59 kg up to 91 kg / 20 - 200 lbs
User age	From 1 to 13 years
Product weight	From 2.7 kg up to 8.6 kg 6 - 19 lbs
Harness	Including
Accessories	Head support Mobile activity tray
Colors	Lilac, Aqua and Gray 

Required user measurement (Total height)	
	Mobile activity tray
Head support	
	Compatible with Special Tomato strollers

## BENEFITS

They are intended for children with mild to moderate positioning needs.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

