

VOYAGE



The ZIPPIE Voyage® early intervention device combines the convenience of a quick-fold, lightweight stroller frame with incredibly versatile seats and positions, providing more options to support growing babies and their caregivers.

Tilt and recline. Improves positioning to help with feeding, digestion, respiratory function and visual orientation. **Image 1.**

Reversible seat with tilt and recline. **Image 2.**

Accessories for daily life. **Image 3.**

The upper part and the frame are separated for practical, easy transport and allows the assistant good handling. **Image 4.**



Image 1



Image 2



Image 3



Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	6" a 13"
Seat depth	5" a 13"
Product weight	13 kg / 28 lbs
Backrest height	7" a 17"
Weight capacity	34 kg / 75 lbs
Recline angle	5° a -30°
Tilt angle	5° Anterior 45° Posterior
Front wheels	7"
Rear wheels	12"
Accessories	Optional
Colors	On base and hoods

Versatile seats
with built-in growth



Options designed
to make everyday
tasks easier



Sistema seguro
para accesorios
de asistencia



BENEFITS

A special needs stroller that provides mild to moderate positioning support

Discover our products **ONE**



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com



ver_2025.11.11E

