

XT4

magic
Mobility



Upgrade to driving the world's most advanced four-wheel drive power wheelchair. The XT4 builds on Magic Mobility's impressive experience in delivering extreme off-road 4WD power wheelchair performance. The XT4 boasts incredible stability, range, comfort, and control over even the most challenging of terrains.

The XT4's patented front and rear suspension. **Image 1.**

High-capacity battery. **Image 2.**

Low-pressure knobby tires. **Image 3.**

Offers three standard colors, three specialized prints, and carbon fiber. **Image 4.**



Image 1



Image 2



Image 3




Image 4



Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	13" up to 23"
Seat depth	13" up to 24"
Product weight	From 320 lbs / 145 kg
Weight capacity	399 lbs / 348 lbs Elevated or tilted
Backrest height	16.5" up to 26"
Maximum range	19 miles
Maximum speed	6.2 mph
Wheels	14"
Turning radius	55.5"
Total width	27.5"



- Seat lift and tilt 12" / 45°
- 95-180° Legrest elevation
- Tilt 50°
- 95-160° Recline

BENEFITS

For users with partial or total reduced mobility of lower limbs.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product
Contact your local Loh Medical representative**

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

