

# Bariatric Wheelchair



A strong steel wheelchair made for patients that require extra weight capacity. It is comfortable, foldable, patients with active lifestyles seeking a wheelchair with high performance. Ideal for: Patients requiring additional space and robust support in a wheelchair for daily use. For users weighing over 150kg.

Foldable frame for easy transportation and storage. **Image 1.**

Desk length folding armrests, PVC armrests. **Image 2.**

Reinforced frame with double cross and stabilizing bar on the back to support greater user weight. **Image 3.**

Aluminum footrest with heel band. **Image 4.**



Image 1



Image 2



Image 3



Image 4




# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Frame	Powder-coated steel folding
Seat width	22" and 24"
Seat depth	18"
Product weight	24 kg
Weight capacity	150 kg
Front wheel	8"x2" solid
Rear wheel	24" solid
Armrest	PVC folding desk
Footrest	Removable and height adjustable
Fabric	Extra strong fire resistant nylon with foam
Medicare Code	K0006 & K0007



### Shipping Measurements

A 33"

B 11"

C 37"

Weight 28 kg.

## BENEFITS

For people who require a light and safe transfer chair.

## Discover our products



### DESIGNED FOR OUR TERRITORIES


Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2025.11.11E

