

# Power Chair with electric tilt & recline

**FELIX<sup>RT</sup>**



**ONE**

The Felix<sup>®</sup> RT is a motorized wheelchair with a rigid frame and electronic tilt and recline system that offers maximum independence and comfort. Designed for users with neurological damage, it allows for effortless postural adjustments. Its injection-molded foam backrest and seat provide continuous support, helping to maintain proper and comfortable posture at all times.

Dynamic controller + lights + Actuator and USB. **Image 1.**

Front and rear lights. **Image 2.**

Electronic tilt from 0° to 20°. **Image 3.**

Electronic recline from 90° to 140°. **Image 4.**



Image 1




Image 2



Image 3




Image 4




# Keeping you in mind. More options.

Quality products with the required accessories at a great price.

## SPECIFICATIONS

Frame	Rigid Blue steel
Seat Width x Depth	18" x 17"
Backrest height	20.5"
Backrest type	Padded Fabric
Product weight	102.4 kg / 226 lbs
Weight capacity	136 kg / 300 lbs
Wheels	13.5" x 2.5" Solid   8.5" x 2.5" Solid
Motors	2 x 340 W, Electromagnetic Brake
Max. Speed	10 kmh
Endurance mileage	45 km
Battery	24 V, 52 Ah Lead-Acid
Suspension	Rear
HCPCS	K0013



**BENEFITS** For users with reduced mobility partial or total lower limbs.

## Discover our products



### DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



### FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

**For more information on obtaining this product  
Contact your local Loh Medical representative**

Find your local representative at [www.lohmedical.com](http://www.lohmedical.com)   

ver\_2026.04.11E

