

Jay Lite Pro

JAY
Superior Clinical Seating.



Highly active self-propellers need seating that keeps up and boosts propulsion – reliable support during every transfer, every push, every challenge. The JAY LitePro cushion meets that demand with uncompromising performance. Using a state-of-the-art multi-zone elastomer insert, the LitePro provides pressure relief, breathability, and stability in a durable, low-profile, ultra lightweight cushion.

At just 2 lbs. (0.9 kg) and only 2" (51 mm), the JAY LitePro is one of the lightest and lowest-profile adjustable skin protection cushions on the market. Its slim design integrates seamlessly into the wheelchair without adding bulk.

Figure 1

Skin health starts with moisture control. With 45% better humidity transfer than a traditional air cushion, the JAY LitePro cushion helps keep skin cooler and drier throughout the day.

Figure 2

Integrated, interchangeable pressure management pods allow personalization of the level of pressure relief to individual client needs. **Figure 3**

The patent-pending elastomer insert allows for up to 1.6" (40 mm) of immersion, helping to redistribute pressure away from high-risk bony areas and support healthier, better-protected skin. **Figure 4**



Image 1



Image 2



Image 3




Image 4

Keeping you in mind. More options.

For those who require cutting-edge technology to keep up with the demands of their active life.

SPECIFICATIONS

Seat width	14" to 20"
Seat depth	14" to 20"
Product height	2"
Product weight	Less than 2 lbs

One of the lightest and lowest-profile cushions on the market



Combines long-lasting performance with simple, low-maintenance cleaning.

The insert maintains consistent skin protection, even during frequent transfers



BENEFITS

Primarily designed for active users of ultralight manual wheelchairs who self-propel and require a lightweight, stable cushion with a high level of skin protection.

Discover our products



DESIGNED FOR OUR TERRITORIES

Our products are designed to meet the needs of each of our territories, offering quality and affordable prices.



FUNCTIONAL VERSALITY

We have a variety of products ranging from a wide range of manual and power wheelchairs and mobility scooter.

For more information on obtaining this product
Contact your local Loh Medical representative

Find your local representative at www.lohmedical.com   

ver_2025.11.11E

